AUDIT® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems

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ENGLISH

INTENDED USE

The Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems is assayed quality control material consisting of two sets of five levels of human and bovine serum albumin. Each level of Set 1 contains the following analytes: CEA, Ferritin, Free PSA, FSH, LH, Progesterone, Prolactin, Total PSA, Total T3, Total T4, TSH. Each level of Set 2 contains the following analytes: Free T4. These five levels of each set demonstrate a linear relationship to each other for their respective analytes¹. It is intended to simulate human patient serum samples for purpose of determining linearity, calibration verification and verification of reportable range for CEA, Ferritin, Free PSA, FSH, LH, Progesterone, Prolactin, Total PSA, Total T3, Total T4, TSH and Free T4.

This product is intended for use with quantitative assays on the indicated analyzer provided in the labeling. The Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems should not be used for calibration or standardization of the CEA, Ferritin, Free PSA, FSH, LH, Progesterone, Prolactin, Total PSA, Total T3, Total T4, TSH, and Free T4 assays. The Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems is for In Vitro Diagnostic use only.

SUMMARY AND PRINCIPLE

As defined in the Clinical Laboratory Improvement Amendments of 1988 (CLIA) by the Centers for Medicare and Medicaid Services (CMS) and the Centers for Disease Control (CDC), each laboratory must validate each test method's AMR at least every six months as well as following changes in lots of analytically critical reagents or major system components². Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems may be used as one would use human whole blood to verify and validate the AMR.

WARNINGS AND PRECAUTIONS

Because this product is of human origin, it should be handled as though capable of transmitting infectious diseases. Each serum, plasma or whole blood donor unit used in the preparation of this material was tested by United States Food and Drug Administration (FDA) approved methods and found to be negative for antibodies to HIV and HCV and nonreactive for HBSAg. Because no test method can offer complete assurance that HIV, hepatitis B virus, and hepatitis C virus or other infectious agents are absent, this material should be handled as though capable of transmitting infectious diseases. This product may also contain other human source material for which there is no approved test. The FDA recommends such samples be handled at the Centers for Disease Control's Biosafety Level 2.

This product contains less than 0.1% sodium azide that may react with lead and copper plumbing to form potentially explosive metal azides. On disposal, flush with a large volume of water to prevent azide build-up.

Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems is intended solely for the purpose of in vitro diagnostic use as described on the label. Audit® MicroControls™, Inc. will not be liable for any unclaimed damages arising from any other usage.

MATERIALS PROVIDED

Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems, 10 x 5 mL

STORAGE AND STABILITY

Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for its assays, Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems has a reconstituted stability of up to 5 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

PROCEDURE

Follow the manufacturer's instructions provided for quality control and for verifying and validating the AMR. Verify that the lot number on each vial matches the package insert. To avoid evaporation, do not leave the vial uncapped. Q.C. requirements should be performed in conformance with local, state and/or federal regulations or accreditation requirements. Calibration verification linearity material should be run³:

- every six (6) months.
- when a complete change of reagents for a procedure is introduced.
- when there is major preventive maintenance or replacement of critical parts that may influence test performance.
- when control materials reflect an unusual trend or shift, or are outside of the laboratory's acceptable limits.
- when the laboratory's established schedule for verifying the reportable range for patient test results requires more frequent calibration verification.

INSTRUCTIONS FOR USE

- Remove a vial from the package, reconstitute only vials B through E with 5ml of deionized water, and gently swirl occasionally for 10 minutes. Do not shake. Do not mix mechanically. Vial A does not need to be reconstituted.
- Refer to instrument or assay instruction manual for verifying and validating the AMR.
- After sampling, replace stopper and return to original package at 2-8°C to obtain the maximum open vial stability.

CALCULATIONS OF RESULTS

Levels B through D in each set of Audit® MicroCV™ Immunoassay Linearity Set for use with Architect® i Immunoassay Systems are prepared in a manner such that an equal distance exists between each consecutive level. This dilution scheme is consistent with the NCCLS recommendations for preparing linearity sets. Level A has been prepared as a "zero level", while Level E has been prepared to challenge the auto dilution function of the analyzer.

U.S. customers only - Once each vial of levels B through D has been tested, raw data may be entered via the AUDITOR QC Program at www.auditmicro.com. An on-line graph showing actual values versus predicted values for each analyte is then available to print, along with slope and intercept data. Call (866) 25-AUDIT for more information.

LIMITATIONS OF THE PROCEDURE

Make sure that each vial is brought to room temperature before testing. If the contents of any of the vials become frozen, discard all vials and request a replacement set, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

EXPECTED VALUES

Each lot of product is manufactured such that a linear relationship exists among levels B - D. The analyte concentrations in this insert were derived from multiple replicate analyses. Actual results obtained may vary depending on instrumentation, methodology and assay temperature. Results may also be dependent on the accuracy of the instrument/reagent system calibration. The degree of acceptable non-linearity is an individual judgment based on methodology, clinical significance and medical decision levels of the test analyte. The material and information presented here in no manner constitutes an overruling of any federal, state or other regulatory body's regulations and/or guidelines.

ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K831M-5	Immunoassay Linearity Set for use with Architect® i Immunoassay Systems	10 x 5 mL

Distributed by AUDIT MicroControls, Inc. - U.S. customers only please call (866) 252-8348 or www.auditmicro.com

¹ Dilution schemes are based on guidelines provided by The National Committee for Clinical Laboratory Standards (NCCLS) in approved guideline EP6-A, "Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline", April 2003.

² Federal Register 42 CFR Part 493, Department of Health and Human Services, January 24, 2003, p.3690.


³ Federal Register 42 CFR Part 493, Department of Health and Human Services, January 24, 2003; s493.1255, (b) (1) (ii).

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