# AUDIT<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control

**REF** K062M-6

6 x 3 mL

**LOT** 06245

 11/30/2012

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## ENGLISH

### INTENDED USE

Audit<sup>®</sup>MicroLQ<sup>™</sup> Ammonia/Ethanol Control is an assayed, stable, ready-to-use liquid, bi-level reference control for use with assays designed to quantitate Ammonia and Ethanol. It is intended to simulate human patient samples for the purpose of monitoring the precision of laboratory testing procedures for Ammonia and Ethanol assays. When used for quality control purposes, it is recommended that each laboratory establish its own means and acceptable ranges and use the values provided only as guides.

This product is intended for use with quantitative assays on the indicated analyzer provided in the labeling. The Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control should not be used for calibration or standardization of the Ammonia and Ethanol assay. The Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control is for In Vitro Diagnostic use only.

### SUMMARY AND PRINCIPLE

Good laboratory practices require that stable reference materials be used to verify the accuracy and precision of testing methods and techniques. Audit<sup>®</sup>MicroLQ<sup>™</sup> Ammonia/Ethanol Control may be used as one would use human serum to obtain the stated Ammonia and Ethanol values. This control will assist in the evaluation of proper performance of Ammonia and Ethanol assays.

### REAGENTS

Audit<sup>®</sup>MicroLQ<sup>™</sup> Ammonia/Ethanol Control is prepared using bovine serum albumin with Ammonia and Ethanol added to achieve the desired concentration levels. Audit<sup>®</sup>MicroLQ<sup>™</sup> Ammonia/Ethanol Control is a ready to use liquid product manufactured according to standard quality control procedures. The manufacturer guarantees stability and consistency of this product.

### WARNINGS AND PRECAUTIONS

Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control is intended solely for in vitro diagnostic use for the purpose described on the labeling. Audit<sup>®</sup> MicroControls<sup>™</sup>, Inc. shall not be liable for any unclaimed damages arising from any other usage.

### MATERIALS PROVIDED

Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control, 6 x 3 mL

### STORAGE AND STABILITY

Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control is stored at 2-8°C and will remain stable in the unopened vial until the expiration date. After opening, the contents should be used according to the instrument manufacturer's instructions and immediately returned to 2-8°C.

When used to monitor the precision of laboratory testing procedures for Ammonia and Ethanol assays, Audit<sup>®</sup> MicroLQ<sup>™</sup> Ammonia/Ethanol Control has an open vial stability of up to 30 days under the proper storage conditions. Leaving the vial uncapped, or prolonging its time at room temperature, will void this open vial stability claim. Make sure the contents of the vial are well mixed before use.

## PROCEDURE

Follow the manufacturer's instructions provided for Ammonia and Ethanol procedures. Verify that the lot number on the vial matches the assay sheet. To avoid evaporation, do not leave the vial uncapped. Controls should be run:

1. daily, in conjunction with patient samples.
2. as recommended by the instrument manufacturer.
3. as required by the relevant regulatory agency.

### INSTRUCTIONS FOR USE

1. Remove a vial from the package and mix by inversion. Do not shake. Do not mix mechanically.
2. Insert the needle of a syringe into the septum of the vial for sample transfer.
3. Refer to instrument or assay instruction manual for analyzing control material.
4. After sampling, replace stopper and return to original package for maximum open vial stability at 2-8°C.

### EXPECTED VALUES

The performance range for each level, based on data by combining estimates of assay variance as determined by participating laboratories using approved FDA instruments and reagents, is provided below. Average values obtained in the laboratory should fall within the performance range although the recovery may not be identical with the mean value listed. Variation between labs will be greater than the precision for any one instrument. Accuracy and precision depend on differences in equipment, reagents, supplies and techniques. Therefore, a lab must establish its own acceptable target values and ranges.

### LIMITATIONS OF THE PROCEDURE

Variations in instruments and in the temperature of the testing material may result in accuracy differences. Make sure that the vial is brought to room temperature before testing. If the liquid in the vial becomes frozen, discard and use another vial, as the results will not be valid.

Dispose of any discarded materials in accordance with the requirements of your local waste management authorities.

### ORDERING INFORMATION

PRODUCT NUMBER	PRODUCT DESCRIPTION	PRODUCT PACKAGING
K062M-6	Ammonia/Ethanol Control	6 x 3 mL

Distributed by AUDIT MicroControls, Inc. - U.S. customers only please call (866) 252-8348 or [www.auditmicro.com](http://www.auditmicro.com)

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